







SYMBOL GLOSSARY

This symbol glossary contains all Symbols used across all media generated by Ostoform Ltd.

REGULATORY

	MEDICAL DEVICE SYMBOL Indicates the item is a medical device.		BATCH CODE Indicates the manufacturer's batch code so that the batch or lot can be identified.
	CATALOGUE NUMBER Indicates the manufacturer's catalogue number so that the medical device can be identified.		UNIQUE DEVICE IDENTIFIER Indicates a carrier that contains unique device identifier information.
	DO NOT RE-USE Indicates a medical device that is intended for one single use only.		INSTRUCTIONS FOR USE Indicates the need for the user to consult the instructions for use.
	USE BY DATE Indicates the date after which the medical device is not to be used.		PATIENT INFORMATION WEBSITE Indicates a website where a patient can obtain additional information on the medical product.
	TEMPERATURE LIMIT Indicates the temperature limits to which the medical device can be safely exposed.		UPPER LIMIT OF TEMPERATURE Indicates the upper limit of temperature to which the medical device can be safely exposed.
	KEEP AWAY FROM SUNLIGHT Indicates a medical device that needs protection from light sources.		KEEP DRY Indicates a medical device that needs to be protected from moisture.
	CAUTION Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.		MANUFACTURER Indicates the medical device manufacturer.

REGULATORY - EU SPECIFIC



CE MARKING

On commercial products, the letters CE mean that the manufacturer or importer affirms the good's conformity with European health, safety, and environmental protection standards

REGULATORY - UK SPECIFIC

UKRP

UK RESPONSIBLE PERSON


Indicates the authorised representative who acts on behalf of a non-UK manufacturer within the Great Britain market (England, Wales, and Scotland).



UK CONFORMITY ASSESSED


This marking is a conformity mark that indicates conformity with the applicable requirements for products sold within Great Britain

OSTOFORM MOLDABLE SEAL WITH FLOWASSIST PROTECTION



OSTOFORM SEAL

Indicates on packaging, labelling, Information for use and marketing material that adjacent information pertains to the Ostoform Seal with FLOWASSIST Protection.



SEAL SIZE

Indicates on packaging, labelling, Information for use and marketing material which size Ostoform Seal with FLOWASSIST Protection the adjacent material/text is referring to. Size will be indicated by text in the center of the Seal, along an arrow.



CORRECT ACTION

Displayed as part of or adjacent to an instructional image. Indication that the displayed image is the correct outcome or action.



INCORRECT ACTION


Displayed as part of or adjacent to an instructional image. Indication that the displayed image is the incorrect outcome or action.



CUT WITH SCISSORS


Indicates in the Information for use and marketing material when/how you should cut to accommodate the Ostoform FLOWASSIST Spout using scissors.

OSTOFORM FLOWASSIST TWO PIECE POUCHING SYSTEM




BASEPLATE

Indicates on packaging, labelling, Information for use and marketing material that adjacent information pertains to the Ostoform FLOWASSIST Flat Baseplate.




DRAINABLE POUCH

Indicates on packaging, labelling, Information for use and marketing material that adjacent information pertains to the Ostoform FLOWASSIST Drainable Pouch.



UROSTOMY POUCH

Indicates on packaging, labelling, Information for use and marketing material that adjacent information pertains to the Ostoform FLOWASSIST Urostomy Pouch.



CLOSED POUCH

Indicates on packaging, labelling, Information for use and marketing material that adjacent information pertains to the Ostoform FLOWASSIST Closed Pouch.



CORRECT ACTION

Displayed as part of or adjacent to an instructional image. Indication that the displayed image is the correct outcome or action.



INCORRECT ACTION

Displayed as part of or adjacent to an instructional image. Indication that the displayed image is the incorrect outcome or action.



TIME INDICATED

Indicates the amount of time a step or action should take to complete. The time will always be specified in writing on the left-hand side of the clock face.




BASEPLATE SIZE

Indicates on packaging, labelling, Information for use and marketing material which size Ostoform FLOWASSIST Baseplate the adjacent material/text is referring to. Size will be indicated by text in the center of the Baseplate, along an arrow.




CUT WITH SCISSORS

Indicates on packaging where you should cut the Ostoform FLOWASSIST Baseplate using scissors.



CLIC RING SIZE

Indicates on packaging, labelling, Information for use and marketing material which size Clic Ring the adjacent material/text is referring to. Size will be indicated by text and colour in the center of the Clic Ring.



FLAT BASEPLATE

Indicates on packaging, labelling, that this product is the Ostoform FLOWASSIST Flat Baseplate.



CONVEX BASEPLATE

Indicates on packaging, labelling, that this product is the Ostoform FLOWASSIST Convex Baseplate.