


Application Examples






Filling Gaps on Uneven Skin Surfaces

Use as a packing agent with either:

- FLOWASSIST Seal
- FLOWASSIST 2-Piece Pouch
- Other pouching systems and accessories





Loop Ileostomies

Use either :

- Ostoform Crest + FLOWASSIST Seal
- 2 x Ostoform Crests





Oval Shaped Stomas

Use either :

- Ostoform Crest + FLOWASSIST Seal
- Ostoform Crest + other accessories





Large / Post-Op Stomas

Use with the FLOWASSIST Seal to make an extra-large seal.

- Ostoform Crest + FLOWASSIST Seal

